



Implant System






INDEX

AISER	
• Company	2
• Biosyn-D	4
• Warranty	6
IMPLANT SYSTEM	8
• Tytan	10
• Themys	12
• Hyperion	14
• Ceos	16
CONICAL LINK SYSTEM	18
• Multi-Unit Abutment	20
• Esthetic Abutment Incisal	22
• Esthetic Abutment Molar	24
• Healing Abutment	26
• Ball Attachment	30
• Flat Attachment	32
• Ti Base Abutment	36
• MUA Ti Base - MUA Healing Cap	38
• Screws	40
IMPRESSION COMPONENTS	42
• Scan Body	44
• Transfer Open Tray	46
• Analogs	48
CLINICAL INSTRUMENTS	
• Surgical Kit	50
• Surgical Kit - Advanced	52
• Surgical Kit - Zygomatic	56
• Prosthetic Kit	58

AISER® Dental Implants are the result of continuous research, rigorous quality control, and advanced clinical validation.

The entire production chain is fully traceable through the UTD (Unique Traceability Document), providing assurance to both patients and professionals worldwide.

AISER® Dental Implants are designed for long-term performance. Backed by a lifetime guarantee.



RESEARCH AND CLINICAL TESTS

AISER® Implant System is the result of in-depth research focused on addressing the real clinical challenges practitioners face every day. Our mission is to deliver technologically advanced solutions that enable clinicians to provide the highest standard of care to their patients.

Our R&D division collaborates closely with clinicians, universities, and leading institutions worldwide.

Our core research domains include:

- Advanced surface engineering of titanium and titanium alloys to enhance biocompatibility, osteoinductivity, osteoconductivity, and mechanical durability
- Plasma sputtering technologies enabling antimicrobial surface activity
- High-performance conical implant–prosthetic connections delivering stable engagement at 30 N torque
- Miniaturized and standardized implant systems designed to ensure consistent, reliable outcomes across all clinical indications — with the simplicity of a single connection platform, independent of diameter or implant line





Biosyn-D

AISER proprietary titanium surface treatment represents a breakthrough in implant surface technology.

It produces a stable, highly oxidized, hydrophilic, nanoporous surface — delivering exceptional biological and mechanical performance.

The implant design, combined with a high-roughness surface profile, ensures optimal primary stability.


Enhanced surface **bioactivity**, characterized by osteoinductive activity, fosters rapid and high-quality osseointegration, ensuring consistent and predictable secondary stability.

BIOSYN-D

Swiss precision in biotechnologies.



Study on the surface performed
at the University of Geneva



5.00 µm

S4800 10.0kV 9.3 mm x10.0k SE (M)

The UTD (Unique Traceability Document) ensures full traceability of every product, covering all production stages and logistics — from manufacturing to delivery.

Our quality control processes are among the most strict in the industry, combining high-precision dimensional verification using laser-based systems, with advanced EDS/SEM chemical and topographical analyses.

This level of traceability and systematic quality assurance is what enables AISER® to provide a lifetime warranty across its product range.


AISER® Implants are engineered for lifetime performance.

AISER offers a **lifetime warranty** across its product range.



Mechanical study performed
at the Polytechnic University of Milan

Reliability you can trust over time.



AISER® Implant System

Reliability and exceptional bioactivity.

AISER® Implant System is engineered not to compromise.

High-grade titanium ensures exceptional **mechanical reliability**.

Enhanced **bioactivity performance** through **AISER® Biosyn-D** oxidative etching treatment.

Reliability and **innovation** — delivering an implant system designed to perform across all clinical indications.

AISER® Implant System: Tytan, Themys, Ceos.

// **Characteristics:**

All components are created with the most **resistant and reliable** biocompatible titanium alloys.

The implant's surface is etched and oxidized through **AISER® Biosyn-D** treatment, resulting in a surface with high grade nano-porosity, nano-roughness, hydrophilicity, osteoinductivity.

Implants are biphasic with an internal hexagonal **conical** coupling and a **Universal Prosthetic Platform** accross all implant lines and diameter for unmatched simplicity.


Implants are designed to achieve primary stability, and to guarantee a practical insertion.



TYTAN



HYPERION



THEMYS



TYTAN

Control · Stress Management · Stability



AISER® Implant System

TYTAN

Tytan implant line is specifically engineered for **D1** and **D2** bone types.

Its **dual-helix** geometry provides enhanced primary stability while optimizing insertion dynamics and reducing friction within the osteotomy.

As with all AISER implants, Tytan is equipped with an **hexagonal conical connection** and a **Universal Prosthetic Platform** across **all diameters**.



Universal Prosthetic Platform
ONE CONNECTION
Across all Ø and implant lines



Ø 3,0 mm	Ø 3,5 mm	Ø 3,9 mm	Ø 4,3 mm	L
TY30W06	TY35V06	TY39GR06	TY43G06	6 mm
TY30W07	TY35V07	TY39GR07	TY43G07	7 mm
TY30W08	TY35V08	TY39GR08	TY43G08	8 mm
TY30W10	TY35V10	TY39GR10	TY43G10	10 mm
TY30W11	TY35V11	TY39GR11	TY43G11	11 mm
TY30W13	TY35V13	TY39GR13	TY43G13	13 mm
TY30W15	TY35V15	TY39GR15	TY43G15	15 mm
TY30W18	TY35V18	TY39GR18	TY43G18	18 mm
TY30W22	TY35V22	TY39GR22	TY43G22	22 mm
TY30W24	TY35V24	TY39GR24	TY43G24	24 mm
TY30W26	TY35V26	TY39GR26	TY43G26	26 mm



Ø 5,0 mm	Ø 5,5 mm
TY50B06	TY55Y06
TY50B07	TY55Y07
TY50B08	TY55Y08
TY50B10	TY55Y10
TY50B11	TY55Y11
TY50B13	TY55Y13
TY50B15	N/A





THEMYS

Adaptability · Performance · Stability



AISER® Implant System

THEMYS

Themys is the most **versatile** implant line within the **AISER** portfolio, engineered to perform **reliably** across **all bone types**.

Its **single-helix** geometry with a **wide pitch** and **aggressive thread** profile enables efficient cutting, reduced bone compression, and **outstanding primary stability**.

As with all AISER implants, Themys is equipped with an **hexagonal conical connection** and a **Universal Prosthetic Platform** across **all diameters**.



Universal Prosthetic
Platform
ONE CONNECTION
Across all Ø
and implant lines



Ø 3,0 mm	Ø 3,5 mm	Ø 3,9 mm	Ø 4,3 mm	L
TH30W06	TH35V06	TH39GR06	TH43G06	6 mm
TH30W07	TH35V07	TH39GR07	TH43G07	7 mm
TH30W08	TH35V08	TH39GR08	TH43G08	8 mm
TH30W10	TH35V10	TH39GR10	TH43G10	10 mm
TH30W11	TH35V11	TH39GR11	TH43G11	11 mm
TH30W13	TH35V13	TH39GR13	TH43G13	13 mm
TH30W15	TH35V15	TH39GR15	TH43G15	15 mm
TH30W18	TH35V18	TH39GR18	TH43G18	18 mm
TH30W22	TH35V22	TH39GR22	TH43G22	22 mm
TH30W24	TH35V24	TH39GR24	TH43G24	24 mm
TH30W26	TH35V26	TH39GR26	TH43G26	26 mm



Ø 5,0 mm	Ø 5,5 mm
TH50B06	TH55Y06
TH50B07	TH55Y07
TH50B08	TH55Y08
TH50B10	TH55Y10
TH50B11	TH55Y11
TH50B13	TH55Y13
TH50B15	N/A





HYPERION

Mechanical Engagement · Bone Compaction · Stress Management



HYPERION

Hyperion represents the most advanced implant line within the AISER portfolio, developed to ensure reliable primary stability, even in demanding clinical conditions.

The **apical thread** design supports immediate **mechanical engagement**, while the progressively evolving thread profile—becoming more **square** towards the coronal region—promotes **controlled bone compaction** and reduces stress on the cortical bone.

As with all AISER implants, Hyperion is equipped with an **hexagonal conical connection** and a **single prosthetic platform** across **all diameters**.



AISER® Implant System



Universal Prosthetic
Platform
ONE CONNECTION
Across all Ø
and implant lines



Ø 3,0 mm	Ø 3,5 mm	Ø 3,9 mm	Ø 4,3 mm	L
TH30W06	TH35V06	TH39GR06	TH43G06	6 mm
TH30W07	TH35V07	TH39GR07	TH43G07	7 mm
TH30W08	TH35V08	TH39GR08	TH43G08	8 mm
TH30W10	TH35V10	TH39GR10	TH43G10	10 mm
TH30W11	TH35V11	TH39GR11	TH43G11	11 mm
TH30W13	TH35V13	TH39GR13	TH43G13	13 mm
TH30W15	TH35V15	TH39GR15	TH43G15	15 mm
TH30W18	TH35V18	TH39GR18	TH43G18	18 mm
TH30W22	TH35V22	TH39GR22	TH43G22	22 mm
TH30W24	TH35V24	TH39GR24	TH43G24	24 mm
TH30W26	TH35V26	TH39GR26	TH43G26	26 mm



Ø 5,0 mm	Ø 5,5 mm
TH50B06	TH55Y06
TH50B07	TH55Y07
TH50B08	TH55Y08
TH50B10	TH55Y10
TH50B11	TH55Y11
TH50B13	TH55Y13
TH50B15	N/A





CEOS

Precision, Stability, Tissue Management.



CEOS

Ceos implant line is designed for extrasinus zygomatic surgery.

The implant apex is treated with **AISER® Biosyn-D** technology, which allows for a **reliable** osseointegration.

As with all AISER implants, Ceos is equipped with an **hexagonal conical connection** and a **Universal Prosthetic Platform** across all diameters.



AISER® Implant System



Universal Prosthetic
Platform
ONE CONNECTION
Across all Ø
and implant lines



Ø 3,5 mm	Ø 4,2 mm	L
CE35V30	CE42G30	30 mm
CE35V32	CE42G32	32,5 mm
CE35V35	CE42G35	35 mm
CE35V37	CE42G37	37,5 mm
CE35V40	CE42G40	40 mm
CE35V42	CE42G42	42,5 mm
CE35V45	CE42G45	45 mm
CE35V47	CE42G47	47,5 mm
CE35V50	CE42G50	50 mm
CE35V52	CE42G52	52,5 mm
CE35V55	CE42G55	55 mm
CE35V57	CE42G57	57,5 mm
CE35V60	CE42G60	60 mm

■ CEOS 3,5:
Ø 3,5 mm Fixture - 4,2 mm Collar

■ CEOS 4,2:
Ø 4,2 mm Fixture - 4,2 mm Collar

Specifications:

AISER® Ceos design: Endosseous, Biphasic, single-coil, with wide pitch and thread

AISER® Conical Link System connection: internal, conical hexagonal connection, with conical and screwed implant-prosthetic junction.

AISER® Biosyn-D surface: nano-porous, nano-rough, highly oxidized, therefore surface osteoinductive, osteoconductive and chemically stable.

AISER® Conical Link System

The reliability of a conical connection and the simplicity of a screwed solution.

AISER® Conical Link System was developed by combining two essential concepts: the **reliability** of a **conical connection** and the ease of use of a **single-connection** system designed for everyday clinical practice.

Practicality and **reliability**, to provide clinicians with only the best.


The extensive range of **AISER® Conical Link System** screws and abutments meets virtually every prosthetic requirement, while preserving the simplicity and efficiency of a **Universal Prosthetic Platform** system.

// **Characteristics**


All components are created with the most **resistant and reliable** of biocompatible titanium alloys.

The Implant-Prosthetic junction is built, for every component, with an internal hexagonal **conical connection of 8°**.

All components have been designed with a **switching platform** concept in mind, ensuring optimal integration with the peri-implant soft tissues.




Multi-Unit
Abutment



Esthetic
Abutment



Healing
Abutment



Multi-Unit Abutment

Conical Connection



MULTI-UNIT ABUTMENT

AISER® Multi Unit Abutments

are designed for single and multiple implant-supported rehabilitation, both with immediate or deferred load.

Their geometry facilitates placement while minimizing soft tissue compression, even in highly angulated restorations.

The characteristic “goblet” profile supports **platform switching**, promoting soft tissue trophism, natural aesthetics, and long-term clinical reliability.

The **conical connection**, combined with a dedicated screw and retrieval mechanism, ensures long-term mechanical stability and ease of placement and removal.

AISER® MUA is available in a wide range of angulations, making it an ideal solution for multiple-unit and full-arch rehabilitations.



AISER® Conical Link System

Ø MUA Platform - 5 mm



0°	
MUA00H1	H1
MUA00H2	H2
MUA00H3	H3
MUA00H4	H4
MUA00H5	H5



60°	
MUA60H1	H1
MUA60H2	H2
MUA60H3	H3
MUA60H4	H4
MUA60H5	H5




17°	30°	45°	53°	
MUA17H1	MUA30H1	MUA45H1	MUA53H1	H1
MUA17H2	MUA30H2	MUA45H2	MUA53H2	H2
MUA17H3	MUA30H3	MUA45H3	MUA53H3	H3
MUA17H4	MUA30H4	MUA45H4	MUA53H4	H4
MUA17H5	MUA30H5	MUA45H5	MUA53H5	H5

Specifications:

The fixture-abutment pairing is an internal hexagonal conical connection with 8° angle.

The conical connection allows stable pairing at **30N** torque on the abutment screw.

We provide all Abutments with their dedicated fixing screw, AISER® Universal Clinic Screw.



Esthetic Abutment Incisal

Conical Connection



AISER® Conical Link System

ESTHETIC ABUTMENT

INCISAL

AISER® Esthetics Abutments are designed for single or multiple implant-prosthetic rehabilitation, with immediate or deferred load.

Their geometry facilitates placement while minimizing soft tissue compression.

The characteristic “goblet” profile supports **platform switching**, promoting soft tissue trophism, natural aesthetics, and long-term clinical reliability.

The **conical connection**, combined with a dedicated screw and retrieval mechanism, ensures long-term mechanical stability and ease of placement and removal.



0°	17°	25°	
AB00H1SI	AB17H1SI	AB25H1SI	H1
AB00H2SI	AB17H2SI	AB25H2SI	H2
AB00H3SI	AB17H3SI	AB25H3SI	H3
AB00H4SI	AB17H4SI	AB25H4SI	H4
AB00H5SI	AB17H5SI	AB25H5SI	H5


Specifications:

The fixture-abutment pairing is an internal hexagonal conical connection with 8° angle.

The conical connection allows stable pairing at **30N** torque on the abutment screw.

We provide all Abutments with their dedicated fixing screw, AISER® Universal Clinic Screw.





Esthetic Abutment Molar

Conical Connection



AISER® Conical Link System

ESTHETIC ABUTMENT


MOLAR

AISER® Esthetics Abutments are designed for single or multiple implant-prosthetic rehabilitation, with immediate or deferred load.

Their geometry facilitates placement while minimizing soft tissue compression.

The characteristic “goblet” profile supports **platform switching**, promoting soft tissue trophism, natural aesthetics, and long-term clinical reliability.

The **conical connection**, combined with a dedicated screw and retrieval mechanism, ensures long-term mechanical stability and ease of placement and removal.



0°	17°	25°	
AB00H1SM	AB17H1SM	AB25H1SM	H1
AB00H2SM	AB17H2SM	AB25H2SM	H2
AB00H3SM	AB17H3SM	AB25H3SM	H3
AB00H4SM	AB17H4SM	AB25H4SM	H4
AB00H5SM	AB17H5SM	AB25H5SM	H5

Specifications:

The fixture-abutment pairing is an internal hexagonal conical connection with 8° angle.

The conical connection allows stable pairing at **30N** torque on the abutment screw.

We provide all Abutments with their dedicated fixing screw, AISER® Universal Clinic Screw.





Healing Abutment



AISER® Conical Link System

HEALING ABUTMENT

AISER® Healing Abutment is designed for the post-surgical healing phase, prior to deferred prosthetic loading.

Their geometry facilitates placement while minimizing soft tissue compression.

The characteristic “goblet” profile supports **platform switching**, promoting soft tissue trophism, natural aesthetics, and long-term clinical reliability.



H1	H2	H3	H3	H3
HAH1	HAH2	HAH3	HAH4	HAH5

AISER® Conical Link System

The reliability of a conical connection and the simplicity of a screwed solution.





Ball
Attachment



Flat
Attachment



Ball Attachment

Conical Connection



AISER® Conical Link System

BALL ATTACHMENT

AISER® Ball Attachment is designed to retain removable prostheses, using a nylon (PTFE) retentive component.


Their geometry facilitates placement while minimizing soft tissue compression.

The characteristic “goblet” profile supports **platform switching**, promoting soft tissue trophism, natural aesthetics, and long-term clinical reliability.

The **conical connection**, ensures long-term mechanical stability and ease of placement and removal.



Ball Attachment	
BAH1	H1
BAH2	H2
BAH3	H3
BAH4	H4
BAH5	H5



Ball Attachment
Ball Attachment and Flat Attachment Kit
192ACN

Specifications:

The fixture-abutment pairing is an internal hexagonal conical connection with 8° angle.

The conical connection allows stable pairing at **30N** torque on the abutment screw.





Flat Attachment

Conical Connection



AISER® Conical Link System

FLAT ATTACHMENT

AISER® Flat Attachment is designed to retain removable prostheses, using a nylon (PTFE) retentive component.


Their geometry facilitates placement while minimizing soft tissue compression.

The characteristic “goblet” profile supports **platform switching**, promoting soft tissue trophism, natural aesthetics, and long-term clinical reliability.

The **conical connection**, ensures long-term mechanical stability and ease of placement and removal.



H1	H2	H3	H4	H5
FAH1	FAH2	FAH3	FAH4	FAH5



Flat Attachment
Ball Attachment and Flat Attachment Kit
192ACN


Specifications:


The fixture-abutment pairing is an internal hexagon with double pairing system: screwed and conical.

The conical connection allows stable pairing at **30Ncm** torque on the fixing screw.


AISER® Conical Link System

The reliability of a conical connection and the simplicity of a screwed solution.






Ti Base
Abutment



MUA
Ti Base



Screws



Ti Base Abutment

Conical Connection



AISER® Conical Link System

TI BASE ABUTMENT

AISER® Ti Base is an adaptable abutment, particularly used in digital workflow and designed for single or multiple rehabilitation, with immediate or deferred load.

Their geometry facilitates placement while minimizing soft tissue compression.

The characteristic “goblet” profile supports **platform switching**, promoting soft tissue trophism, natural aesthetics, and long-term clinical reliability.

The **conical connection**, ensures long-term mechanical stability and ease of placement and removal.



ROTATING

Ø 4,5 mm	
TIBH1R	H1
TIBH2R	H2
TIBH3R	H3
TIBH3R	H4
TIBH3R	H5




NOT ROTATING


Ø 4,5 mm	
TIBH1NR	H1
TIBH2NR	H2
TIBH3NR	H3
TIBH4NR	H4
TIBH5NR	H5

Specifications:

The fixture-abutment pairing is an internal hexagon with double pairing system: screwed and conical.

The conical connection allows stable pairing at **30Ncm** torque on the fixing screw.






MUA Ti Base



MUA TI BASE - MUA HEALING CAP

AISER® Multi-Unit Abutment Ti Base, Ti Base Short and Healing Cap are components meant to be used with **AISER® Multi-Unit Abutment**

They connect to MUA through the **AISER® MUA Ti Base Screw**, excepting for the healing cap which come with a built in screw.



MUA: 00°/17°/30°/45°/53°/60°

MUA
Ti Base

MUA-TIB



MUA: 00°/17°/30°/45°/53°/60°

MUA
Ti Base Short

MUA-TIBSH



MUA
Ti Base Screw

MUA-TIBS



MUA: 00°/17°/30°/45°/53°/60°

MUA Healing
Cap 5 mm

MUA-HC5

MUA Healing
Cap 6 mm

MUA-HC6



Screws



AISER® Conical Link System

SCREWS

- **AISER® Cover Screw**
included with each AISER Implant
- **AISER® Universal Clinical Screw**
Implant abutment screw, included with each abutment
- **AISER® MUA Ti Base Screw**
MUA screw, included with each MUA Ti-base
- **AISER® Abutment Extractor**
Universal extractor for AISER® Conical Link System abutment



Cover Screw	Universal Clinical Screw	Extractor Screw
CS	UCS	EXS



MUA Ti Base Screw MUA-TIBS




AISER® Impression Components

Conometry reliability and simplicity in a screwed solution.






Scan
Body



Transfer
Open Tray



Anologs



Scan Body



AISER® Impression Components

SCAN BODY

AISER® Impression Components

Scan Body line consists of:

- **AISER® Scan Body** is designed for the digital acquisition of implant position, both intraorally and in the laboratory.
- **AISER® Scan Body MUA** is designed for the digital acquisition of the Multi Unit Abutment position, both intraorally and in the laboratory.




Scan Body
Implant

I-SB



Scan Body
MUA

MUA-SB



Transfer Open/closed Tray



AISER® Impression Components

TRANSFER OPEN TRAY

AISER® Open and Closed Tray Transfers are designed for the accurate transfer of implant position and features in conventional impression techniques.




Transfer Open Tray	H1
TOT-SH1	TOT-H1



Transfer Open Tray	H2
TOT-SH2	TOT-H2



Transfer Closed Tray	
TCT-S	TCT



Analogs



AISER® Impression Components

ANALOGS

AISER® Implant Analog is designed to replicate the implant position and features during model fabrication.

AISER® MUA Analog is designed to replicate the Multi Unit Abutment position and features during model fabrication.



Implant Analog

I-AN



Multi-Unit Abutment Analog

MUA-AN





Surgical Kit

AISER® Surgical Kit is designed to provide clinicians with all the instruments required for the surgical placement of **AISER® Implant System** fixtures and the installation of **AISER® Conical Link System** prosthetic components.

The casket is a state-of-the-art tool, which allows to sterilize the surgical kit without the removal of any of the heat-sealed rubber inside.







Surgical Kit - Advanced

AISER® Surgical Kit - Advanced is designed to provide clinicians with all the instruments required for the surgical placement of AISER® Implant System units and the installation of AISER® Conical Link System prosthetic components.

Featuring tap drills, bone mills, and three line of drills, the Advanced Surgical Kit supports clinicians across a wide range of clinical indications. a state-of-the-art tool, which allows to sterilize the surgical kit without the removal of any of the heat







Surgical Kit - Zygomatic

AISER® Surgical Kit - Zygomatic is designed to provide clinicians with all the instruments required for the surgical placement of **AISER® Ceos** zygomatic implants and the installation of **AISER® Conical Link System** prosthetic components.







Prosthetic Kit

AISER® Prosthetic Kit is designed to provide clinicians with all the instruments required for the prosthetic phase of **AISER® Implant System** restorations using **AISER® Conical Link System** components.





AISER® Clinical Instruments



SCREWDRIVERS			
Screwdriver Insert Short	Screwdriver Insert Long	ScrewDriver - Ratchet Short	ScrewDriver - Ratchet Long
SD-IS	SD-IL	SD-RS	SD-RL

MUA 0° DRIVER
MUA 0° Driver Insert
MUA00D-I



IMPLANT DRIVERS - HEX DRIVERS	
Hex Driver - Ratchet Short	Hex Driver - Ratchet Long
HD-RS	HD-RL

IMPLANT DRIVERS - HOLDER DRIVERS	
Cricket Driver Insert Short	Cricket Driver Insert Long
CD-IS	CD-IL

ZYGOMATIC IMPLANT - HOLDER / HEX DRIVER	
Cricket Driver Insert E63 Zygomatic	Hex Driver Insert E63 Zygomatic
CD-IE63	HD-IE63

AISER® Clinical Instruments



Twist Drills 6/15 Short								
Twist Drill Ø2.0 - 6/15	Twist Drill Ø2.4-2.8 - 6/15	Twist Drill Ø2.8-3.2 - 6/15	Twist Drill Ø3.0-3.4 - 6/15	Twist Drill Ø3.2-3.6 - 6/15	Twist Drill Ø3.4-3.9 - 6/15	Twist Drill Ø3.8-4.2 - 6/15	Twist Drill Ø4.2-4.6 - 6/15	Twist Drill Ø4.2-5.0 - 6/15
TD20-615	TD24-615	TD28-615	TD30-615	TD32-615	TD34-615	TD38-615	TD42-615	TD4250-615




Twist Drills 10/18 Mid						
Twist Drill Ø2.0 - 10/18	Twist Drill Ø2.4-2.8 - 10/18	Twist Drill Ø2.8-3.2 - 10/18	Twist Drill Ø3.0-3.4 - 10/18	Twist Drill Ø3.2-3.6 - 10/18	Twist Drill Ø3.4/3.9 - 10/18	Twist Drill Ø3.8-4.2 - 10/18
TD20-1018	TD24-1018	TD28-1018	TD30-1018	TD32-1018	TD34-1018	TD38-1018

AISER® Clinical Instruments



Twist Drills 18/25 Long						
Twist Drill Ø2.0 - 18/25	Twist Drill Ø2.4-Ø2.8 - 18/25	Twist Drill Ø2.8-Ø3.2 - 18/25	Twist Drill Ø3.0-Ø3.4 - 18/25	Twist Drill Ø3.2-Ø3.7 - 18/25	Twist Drill Ø3.4-Ø3.9 - 18/25	Twist Drill Ø3.8-Ø4.2 - 18/25
TD20-1825	TD24-1825	TD28-1825	TD30-1825	TD32-1825	TD34-1825	TD38-1825



Countersink	Tap Drills				Lance Drill
Countersink Drill	Tap Drill 3,55	Tap Drill 4,3	Tap Drill 5,0	Tap Drill 5,5	Lance Drill
CSD	TDS-35	TDS-43	TDS-50	TDS-55	LD

AISER® Clinical Instruments



11 - DRILLS STOP - Drills 2-2.8

Drills Stop 15-22 - Ø2-2.8	Drills Stop 13-20 - Ø2-2.8	Drills Stop 11-18 - Ø2-2.8	Drills Stop 10 - Ø2-2.8	Drills Stop 08 - Ø2-2.8
DS15-28	DS13-28	DS11-28	DS10-28	DS08-28



Universal Extender Insert

UE-I



Multi-Unit Abutment Handle

MUA-HL



Direction Indicator

DIRI



11 - DRILLS STOP - Drills 2.8-4.2

Drills Stop 15-22 - Ø2.8-4.2	Drills Stop 13-20 - Ø2.8-4.2	Drills Stop 11-18 - Ø2.8-4.2	Drills Stop 10 - Ø2.8-4.2	Drills Stop 08 - Ø2.8-4.2
DS15-42	DS13-42	DS11-42	DS10-42	DS08-42



BONE MILL

Bone Mill

BM



Bone Mill Guide

BMG



SECOND LAYER

Universal Ratchet
Adapter

URA

AISER® Clinical Instruments



Zygomatic Drills L. 40		
Twist Drill Ø2.0 - Ø2.5 - Ø3.0	Twist Drill Ø2.8 - Ø3.0	Twist Drill Ø2.8 - Ø3.2
TD20-40Z	TD28-40Z	TD32-40Z

Zygomatic Drills L. 50			Special Drills
Twist Drill Ø2.0 - Ø2.5 - Ø3.0	Twist Drill Ø2.8 - Ø3.0	Twist Drill Ø2.8 - Ø3.2	Ball Burr Ø4.2
TD20-50Z	TD28-50Z	TD32-50Z	BB



Zygomatic Drills L. 60		
Twist Drill Ø2.0 - Ø2.5 - Ø3.0 - L. 60	Twist Drill Ø2.8 - Ø3.0 - L. 60	Twist Drill Ø2.8 - Ø3.2 - L. 60
TD20-60Z	TD28-60Z	TD32-60Z

Zygomatic Drills L. 80			Special Drills
Twist Drill Ø2.0 - Ø2.5 - Ø3.0 - L. 80	Twist Drill Ø2.8 - Ø3.0 - L. 80	Twist Drill Ø2.8 - Ø3.2 - L. 80	Knurled Drill Ø4.2
TD20-80Z	TD28-80Z	TD32-80Z	DKS-Z

AISER® Clinical Instruments



SECOND LAYER

Torque Ratchet


TRT



SECOND LAYER

Millimetric Surgical
Probe

MSP



SCREWDRIVER

E63 ScrewDriver


SD-E63




FIX RATCHET

Millimetric Surgical Probe - Zygomatic

MSP-Z



Swiss Quality at its Peak.





Implant System

Rue des Bains, 35
CH-1205 Genève
Switzerland

geneve@aiserimplants.com

+41 22 556 47 37

aiserimplants.com

